## The Hely & Weber Global Knee\*

ADDRESSES PAIN, IMPROVES LIVES

Surgery is often delayed for Bi & Tricompartmental OA patients:

Conservative care guidelines, patient logistical needs,
necessary medical and behavioral improvements
and the emotional readiness of the patient often cause delays.

The Global OA addresses patient pain, provides stability and assists in mobility, helping patients through these times.



Hybrid Wrap Around Format Closed calf section with a wrap around thigh section allows patients to adjust the size of the thigh to accommodate changes in girth.

### Soft, Comfortable Condyle Pads

Soft, oversized Comfort Foam pads accommodate the inflamed bulbous knee.

#### Thin Removable Trakaderm™ "U" Buttress

This is a "positional buttress" which can be adjusted for each individual patient's need. In addition it can be removed if not required as a therapeutic modality.

### - Patella Sling

The "Sling" allows for directional control of the patella and offers a slight tilt adjustment.

Our "Patella Sling" controls the patella with an inferior to superior directional pull.



## Introducing our Dynamic Patella Sling

By controlling the patella in a more natural position in the superior femoral groove, patients can now manage end range gait and are able to propel up and down stairs with less pain and more stability.



## Soft Durable Lycra Pre-Formed Popliteal Panel

A thin, soft fabric in the popliteal area provides more comfort and compression. This results in less skin irritation, an overall better fit, edema control and higher patient compliance.

View from inside of brace reveals extraordinary craftsmanship and attention to detail.

## **Functional knee stability** patients haven't seen in years!



hinge pockets durability and double as a built in pull up assist.

breathable lycra urethane laminate and spacer fabric.

The ROM hinges are slightly offset to accommodate inflamed bulbous knees.



**Bending and** adjusting the struts allows hinges to better fit misshappen knees and apply slight varus or valgus force to the knee.

**Customizable Range of Motion** 



## SIZING CHART Global Knee™ PDAC Approved L1832/33. L2795 add on code is also recon

SIZE	CALF	Thigh
XS	11.0" - 13.0"	15.0" - 17.0"
S	12.5"- 13.5"	16.0" - 18.0"
M	13.5" - 15.5"	17.5" - 19.5"
L	14.5" - 16.5"	18.5" - 20.5"
XL	16.5" - 19.0"	19.5" - 22.0"
XXL	18.0" - 20.5"	20.5" - 23.0"
XXXL	19.0" - 22.0"	24.0" - 27.0"

Larger and custom sizes for 7640 models available upon request. Calf measured at widest part of calf. Thigh measured 6 inches above patella.



See our application video

orders: 817-453-9767 imsexpertsintake@att.net www.imsexperts.net



All content © 2016 Hely & Weber

# The Hely & Weber GLOBAL KNEE

